

The first Simultaneous Single-tube

NAT solution for 6 Sexually Transmitted Disease pathogens

SET-I

Trichomonas vaginalis
Mycoplasma hominis
Ureaplasma urealyticum
Chlamydia trachomatis
Mycoplasma genitalium
Neisseria gonorrhoeae

SET-II

Treponema pallidum
Haemophilus ducreyi
HSV 1 & 2
Candida albicans

TV

CT

UU

NG

MG

MH


Seegene's DSO™ Technology

STD Detection kit

Just with one tube,
You can get results with an excellent accuracy!



Seegene



Today's premiere NAT solution

Simultaneous analysis of 6 pathogens at once

Infinitesimal amount of pathogens can be detected using PCR method.

PCR yields a higher sensitivity than conventional culture methods or IFA. Infinitesimal amount of pathogens, therefore, may be analyzed.

6 different pathogens can be screened simultaneously using multiplex system.

Screening of 6 different viruses is possible by overcoming limitations of general PCR method.

Seegene's DSO™ technology is employed.

No false positives is created due to the superior features of Seegene's DSO™ primers.

4 additional features are strengthened.

Internal control, Contamination prevention system, positive control size marker, and positive control DNA are of focused features to enhance the quality of the kit.

Multiplex-PCR

1) Cost-effective

Consumables and equipments such as tubes, tips, PCR machines, agarose, electrophoresis. Bigger earnings due to all-in-one test on several pathogens

2) Efficient

Carryover contamination problem arising from repetitive experiments is resolved by reduced trials.

Detection of STD pathogens in clinical samples (Set-I).



M : 100bp ladder / 1 ~ 6 : Patient samples
N : Negative control

- ◀ Internal Control
- ◀ T. vaginalis; 580 bp
- ◀ M. hominis; 502 bp
- ◀ U. urealyticum; 435 bp
- ◀ C. trachomatis; 348 bp
- ◀ M. genitalium; 253 bp
- ◀ N. gonorrhoeae; 214 bp

| Sample No. | Pathogen |
|------------|--|
| 1 | N. gonorrhoeae |
| 2 | M. hominis/ U. urealyticum/ C. trachomatis |
| 3 | M. hominis/ C. trachomatis |
| 4 | U. urealyticum/ C. trachomatis |
| 5 | M. hominis/M. genitalium |
| 6 | T. vaginalis/M. hominis |



.... excellence and innovation in molecular biology