


CENOS

.....engineered prescient

AT PATIENT TESTING (APT)[®] PLATFORM

Rapid molecular testing platform for High Consequence Infectious Disease surveillance and detection



Direct from sample
Whole blood · Swabs · Saliva

30 min
Sample to result

No cold chain
Lyophilised reagents

SYSTEM OVERVIEW: Four components. One complete solution.

CENOS Analyser

Portable 2.76 kg automated RT-qPCR with onboard display results.

XF Reagent System

Direct testing from crude samples. No extraction or cold-chain required.


THEIA Software

Advanced data analysis with amplification curves, Ct values, and audit logs.

Assay Panels

Multiplex detection for human, animal, and zoonotic targets.


HOW IT WORKS: Test to result in <30 minutes - 3 simple steps



Prepare XF reagents



Insert sample and start test



Automated detection and reporting

GPS module

Track and map outbreaks in real-time as they evolve – GPS tagged results, wherever you test.

Dual-channel optical engine

High-stability optical system with bicolour fluorescence detection.

Peltier thermal cycler

Rapid and precise thermal cycling – reliable in any environment

Portable

Lightweight: 2.76kg
Compact: 236mm (L) x 173mm (W) x 197mm (D)



Touch screen interface

Built-in display using a simplified "traffic light" output for non-expert users to see results instantly for immediate decision making.

Carbon-polymer reaction vessel

High surface area-to-volume ratio designed to maximise thermal transfer and PCR performance.

Barcode scanner

Integrated barcode scanner identifies the test, initiates the test run and assigns the results to a unique sample ID ensuring complete traceability.

Sample loading port

Designed for safe and efficient sample handling in any environment, supporting a rapid and safe workflow.

Onboard processing system

Autonomous unit that manages temperature control, signal acquisition, automated signal processing and result generation.

Power-flexible electronics

Runs on mains, vehicle batteries, solar, or portable power banks – truly off-grid capable.

TECHNICAL SPECIFICATIONS



Physical characteristics		Power Input (external DC power source required)	
Device dimensions	240 × 175 × 198 mm (H × W × D)	Nominal input voltage: 15 V DC Operating input voltage range: 14 V to 16.5 V DC Maximum permitted DC input voltage: 16.5 V DC Maximum continuous input current: 12 A Maximum input power: 180 W	
Device weight	2.76 kg		
Handle	Integrated flip up handle for safe and secure carrying	Optical system specifications	
Test Storage	≤ 40 test results	Spectrometer based detection	Blue excitation: 488nm. Collection from 530nm to 580nm Red excitation: 620nm. Collection from 650nm-800nm
Operating conditions		Multiplex	Up to 6 targets
Operating temperature	10 to 35 °C	Excitation sources	LED
Operating humidity	10 to 90%	Software & data handling	
Placement requirements	Even, flat, stable surface; kept away from water, corrosive environments, dust, chemical vapours, and direct sunlight	Software interface	THEIA software
Storage & transport conditions		Operating modes	Online and offline (standalone)
Storage temperature	-10 to 40 °C	Data download options	USB, Ethernet, Wi-Fi
Transport method	Must be transported in original Analyser protective case	Device download format	JSON (JavaScript Object Notation)
Shock & vibration protection	IP67 case includes EVA precision cut foam for protection and stabilisation	THEIA export format	JSON (JavaScript Object Notation) or PDF
Case protection	Waterproof and dustproof in the case (IP67), automatic pressure release valve for transportation	Interface & Connectivity	
Analytical / performance specifications		Interface options	Touch screen display, barcode scanner, USB
Measurement principle	RT-qPCR	Connectivity options	Ethernet, Wi-Fi, GPS
Sample preparation	Crude sample inputs	Accessories & compatible components	
Time to Result	≤ 30minutes	Reaction vessel	Analyser reaction vessel
Reaction vessel compatibility	Only with CENOS Analyser reaction vessel	Key	Analyser key
Thermal ramp rate	>8°C/sec	Transport case	IP67 certified, shock absorbing
Instrument warmup requirements	For temperatures below 10°C allow acclimatisation to temperature to minimally 10°C	Additional accessories	Power supply, ethernet data cable, 90° USB adapter
		Regulatory & compliance	
		Manufacturing standard	Designed and manufactured under ISO 13485 certified QMS (BG Research Ltd)
		Safety & labelling symbols	Refer to symbols and glossary section in instruction manual