

Biochemicals

Catalogue Number	Description	Qty
200-0120	IPTG (Dioxan Free)	1g
200-0120/5	IPTG (Dioxan Free)	5g
200-0190	X-Gal	1g
200-0190/5	X-Gal	5g

SPECIFICATION

Biochemical	Empirical Formula	Molecular Weight	Alternative Name
IPTG (Dioxan Free)	C ₉ H ₁₈ O ₅ S	238.3	Isopropyl-beta-D-thiogalactopyranoside
X-Gal	C ₁₄ H ₁₅ BrClNO ₆	408.6	5-bromo-4-chloro-3-indolyl-beta-D-galactopyranoside



BioGene Limited – 6 The Business Centre, Harvard Way,
Kimbolton, Cambridgeshire, PE28 0NJ, United Kingdom.

+44 (0)3333 222678

biogene@biogene.com

Registered No. 2800327



+44 (0)3333 222567

www.biogene.com

VAT No. GB576854584



5282-QMS-001
ISO 9001



0026