

## BioSizer-1kb™ DNA Ladder

Catalogue Number	Description	Qty
BG-M1000	BioSizer-1kb DNA Ladder	50ug


Bio1Kb DNA Ladder contains 11 linear double strand DNA bands of 500, 1000, 1500, 2000, 2500, 3000, 4000, 5000, 6000, 8000 and 10,000 bases pairs. The intensity of 3000bp band is stronger than the others, which is useful as a reference for band identification.


### ANALYSIS


The marker is evaluated for resolution, intensity, and background on 0.8% agarose gel in 0.5X TBE buffer.


### STORAGE CONDITION


-20°C, stable for one year



 BioGene Limited – 6 The Business Centre, Harvard Way,  
Kimbolton, Cambridgeshire, PE28 0NJ, United Kingdom.

 +44 (0)3333 222678

 biogene@biogene.com

 +44 (0)3333 222567

 [www.biogene.com](http://www.biogene.com)

Registered No. 2800327

VAT No. GB576854584



5282-QMS-001  
ISO 9001



0026